



Technical Information

AUTOSOL® Cutlery Rust Remover

AUTOSOL® Cutlery Rust Remover for effective removal of rust stains from stainless steel cutlery. Activated citric acid resolves stains and discolouration. Especially selected polishing particles remove the solved stains and renew the brilliant gloss of the surface.

Applications

Shake bottle well before use. Apply a small amount of the product to the cutlery and let take effect for a short time. Polish the cutlery with a moist cloth or sponge. Rinse thoroughly with water.

Note

Do not use on gold or silver plated cutlery. Avoid direct food contact. On brushed or matt surfaces use only little pressure.

Technical Data

Form: liquid
Colour: white
Odour: characteristic

Health and safety

Safety Data Sheet is available separately.

Additional Information

Art.-No.: 11 040700

bottle 250 ml

PU: 12



The information provided in this Technical Information is based on the current status of knowledge and experience and is intended to describe the products with respect to possible applications and safety requirements. This information does not represent a legally binding assurance of specific characteristics. A legally binding assurance of specific characteristics or suitability for the concrete purpose of use cannot be derived from our information. Possible trademarks and existing legislation and stipulations are to be observed by the recipient of our products under his/her own responsibility.